

EliteVue

Fibre Optic Macro Otoscope

The EliteVue otoscope is setting the standards.

The EliteVue offers twice the field of view and 66% greater magnification compared to standard otoscopes*.

The EliteVue offers a significantly improved view of the auditory canal and tympanic membrane, enabling more efficient diagnosis and enhanced patient outcomes. Suitable for throat illumination.

- High-performance LED illumination (4,000 Kelvin, CRI = 92)
- More than 5-fold magnification
- Twice the field of view
- Complete view of the tympanic membrane
- Reflection-optimized high performance glass optics
- Wheel to adjust focal depth
- Easy to operate SpecEjec specula ejection device at the rear side of the otoscope to avoid contamination
- Suitable for pneumatic otoscopy (supplied without ball)
- A wide range of power supply options: as part of a diagnostic station, battery operated or with a rechargeable handle.

EliteVue is compatible with ri-scope® L handles and power sources as well as the ri-former wall diagnostic station.



Technical Data*

	EliteVue	ri-scope® L
Field of view, Ø	11 mm (0.59 in)	5,5 mm (0.59 in)
Magnification	5x	3x

* at a distance of 15 mm (0.59 in) between landmark and distal tip of the specula.



Order information – EliteVue Macro Otoscope

Heads

Art. No.: 10510	EliteVue head only, LED, 2.5 V
Art. No.: 10511	EliteVue head only, XL, 2.5 V
Art. No.: 10512	EliteVue head only, LED, 3.5 V
Art. No.: 10512-301	EliteVue head only, LED, 3.5 V, with anti-theft (for ri-former wall diagnostic station)
Art. No.: 10513	EliteVue head only, XL, 3.5 V
Art. No.: 10513-301	EliteVue head only, XL, 3.5 V, with anti-theft (for ri-former wall diagnostic station)

Otoscope Sets

Art. No.: 2200-204	EliteVue otoscope set, LED, 2.5 V, with C-handle for 2 alkaline batteries
Art. No.: 2200-202	EliteVue otoscope set, XL, 2.5 V, with C-handle for 2 alkaline batteries
Art. No.: 2200-203	EliteVue otoscope set, LED, 3.5 V, with C-handle for 1 rechargeable Li-ion battery
Art. No.: 2200-201	EliteVue otoscope set, XL, 3.5 V, with C-handle for 1 rechargeable Li-ion battery

Oto-/Ophthalmoscope Sets

Art. No.: 2210-204	EliteVue oto-/ophthalmoscope L2 set, LED, 2.5 V, with C-handle for 2 alkaline batteries
Art. No.: 2210-202	EliteVue oto-/ophthalmoscope L2 set, XL, 2.5 V, with C-handle for 2 alkaline batteries
Art. No.: 2210-203	EliteVue oto-/ophthalmoscope L2 set, LED, 3.5 V, with C-handle for 1 rechargeable Li-ion battery
Art. No.: 2210-201	EliteVue oto-/ophthalmoscope L2 set, XL, 3.5 V, with C-handle for 1 rechargeable Li-ion battery

Otoscope Sets rechargeable*

Art. No.: 2211-203	EliteVue otoscope set, LED, 3.5 V, with 1 C-handle, 1 rechargeable Li-ion battery and plug-in charger
--------------------	--

Oto-/Ophthalmoscope Sets rechargeable*

Art. No.: 2212-203	EliteVue oto-/ophthalmoscope set, LED, 3.5 V, with 2 C-handles, 2 rechargeable Li-ion batteries and desk charger
--------------------	---

Charger*

Art. No.: 10707	Plug-in charger FOR Li-ion battery 3.5 V ri-accu [®] L for C-handle
Art. No.: 10708	Plug-in charger WITH Li-ion battery 3.5 V ri-accu [®] L for C-handle
Art. No.: 10705	Desk charger ri-charger [®] L 3.5 V/230 V with wall support
Art. No.: 10706	Desk charger ri-charger [®] L 3.5 V/120 V with wall support

*country-specific plugs on request

